

Opilex®

OPILEX® WOUND DRESSING™

OPILEX® WOUND DRESSING is a disposable non-woven adhesive dressing for wound protection. OPILEX® WOUND DRESSING has an easy peeling technology based on hypoallergenic adhesives, and also is designed to prevent secondary infection at the application regions with dual safety technology. OPILEX® WOUND DRESSING is consist of high air permeable pad that rapidly absorbs exudate fluid, and provides a stable tissue healing environment. It uses a soft non-woven fabric for protection of sensitive skin and air permeability.



How to use WOUND DRESSING




Description

Name	OPILEX® WOUND DESSING
Classification	Quasi-drug
Purpose of use	Wound protect
Instruction for use	Non-woven dressing at the application site (PRN)
Ingredient	Non-woven adhesive tape/Polyethylene mesh embo pad
Manufacture	OPILEX Co., LTD.
Country of manufacture	Republic of Korea

Standard specification

Product No.	Size	Packing unit	CTN
OPW-4555	4.5 cm × 5.5 cm	50 EA	2,000 EA
OPW-6070	6 cm × 7 cm	50 EA	2,000 EA
OPW-0910	9 cm × 10 cm	50 EA	2,000 EA
OPW-0915	9 cm × 15 cm	50 EA	2,000 EA
OPW-0920	9 cm × 20 cm	50 EA	1,500 EA
OPW-0925	9 cm × 25 cm	50 EA	1,000 EA
OPW-0930	9 cm × 30 cm	50 EA	1,000 EA
OPW-0935	9 cm × 35 cm	50 EA	1,000 EA



Customer service

Address : 4-8, Gadong 2-gil, Jeonggwan-eup, Gijang-gun, Busan, Republic of Korea
Phone : +82-1577-5264 / E-mail : opilex@opilex.co.kr
Web : www.opilex.co.kr

Contact for exclusive distributor

Phone : +82-51-728-8202 / Fax : +82-51-727-8203
E-mail : opilex@opilex.co.kr





OPILEX CO., LTD.

Opilex®

OPILEX® WOUND DRESSING™

Dressing innovative technologies of OPILEX Co., LTD.

Dressing accuracy of treatment
Dressing adhesive does not stick to hand
Dual removal technology of release paper
Secondary infection protective technology



OPILEX CO., LTD.