



Opilex®

## OPI® FIXING ROLL™

OPI® FIXING ROLL is a disposable non-woven adhesive tape that used for the purpose of fixing gauze, bandage, and catheter. You can use it by cutting it to an appropriate size depending on application regions.

OPI® FIXING ROLL has an easy peeling technology based on a hypoallergenic adhesive, and also is designed to prevent secondary infection at the application regions with dual safety technology.

OPI® Fixing Roll uses a soft non-woven fabric for protection of sensitive skin and air permeability.



### How to use FIXING ROLL




#### Description

Name	OPI® FIXING ROLL
Classification	Quasi-drug
Purpose of use	The fixation of gauze, bandage, and catheter
Instruction for use	Cutting depending on application regions
Ingredient	Non-woven adhesive tape
Manufacture	OPILEX Co., LTD.
Country of manufacture	Republic of Korea

#### Standard specification

Product No.	Size	Packing unit	CTN
OPIFR-02510	2.5 cm × 10 M	1 EA	200EA
OPIFR-0510	5 cm × 10 M	1 EA	200EA
OPIFR-1010	10 cm × 10 M	1 EA	100EA
OPIFR-1510	15 cm × 10 M	1 EA	100EA
OPIFR-2010	20 cm × 10 M	1 EA	60EA
OPIFR-2510	25 cm × 10 M	1 EA	60EA
OPIFR-3010	30 cm × 10 M	1 EA	40EA



#### Customer service

Address : 4-8, Gadong 2-gil, Jeonggwan-eup, Gijang-gun, Busan, Republic of Korea

Phone : +82-1577-5264 / E-mail : opilex@opilex.co.kr

Web : www.opilex.co.kr

#### Contact for exclusive distributor

Phone : +82-51-728-8202 / Fax : +82-51-727-8203

E-mail : opilex@opilex.co.kr



## OPILEX CO., LTD.



# OPI<sup>®</sup> FIXING ROLL<sup>™</sup>

Dressing innovative technologies of OPILEX Co., LTD.

Dressing accuracy of treatment  
Dressing adhesive does not stick to hand  
Dual removal technology of release paper  
Secondary infection protective technology



OPILEX CO., LTD.