


Opilex®

OPILEX® FLEX™

오피렉스 플렉스는 황갈색 폴리우레탄필름 탄력원단을 사용하여 주사 처치 후 상처부위의 압박, 고정을 통한 지혈에 도움을 줍니다. 고탄력 원단의 적용으로 관절부위 또는 움직임이 많은 부위의 사용시에도 신축성에 의한 압박이 유효합니다. 적용되는 지혈 부위의 출혈량에 따라 일반형과 대형을 구분하여 사용할 수 있습니다. 저 자극 점착제를 사용하여 피부트러블 걱정을 줄였으며 30분 이내에 제거하여야 합니다.



How to use FLEX




Description

제품명	오피렉스 플렉스
제품구분	의약외품
효능 효과	상처부위, 환부 등의 보호
용법 용량	환부에 1일 1회 내지 수회 1매씩 갈아붙인다.
전성분	황갈색폴리우레탄점착부직포, 패드(부직포/폴리에틸렌망)
제조원	주식회사 오피렉스/부산광역시 북구 효열로 111 부산지식산업센터 109, 110호
원산지	MADE IN KOREA

Standard specification

Product No.	Size	Packing unit	CTN	EDI급여
OPF-0275	2cm × 7.5cm	50 EA	3,000 EA	K9206003
OPF-0475	4cm × 7.5cm	50 EA	1,500 EA	K9206003
OPF-0510	5cm × 10cm	50 EA	1,000 EA	K9206003



고객문의

주소 : 부산광역시 기장군 정관읍 가동2길 4-8(1F)

Tel : 051-728-8202 / Fax : 051-727-8203

E-mail : opilex@opilex.com / Web : www.opilex.co.kr

Customer service

Address : 4-8, Gadong 2-gil, Jeonggwan-eup, Gijang-gun, Busan, Republic of Korea

Phone : +82-1577-5264 / E-mail : opilex@opilex.co.kr / Web : www.opilex.co.kr

OPILEX CO., LTD.


Opilex® 오피렉스® Opi® 오피아이® 는 주식회사 오피렉스의 의료용품 브랜드입니다.
Opilex® 오피렉스® Opi® 오피아이® is the medical product brand of Opilex.



OPILEX[®] FLEX[™]

Dressing innovative technologies of OPILEX Co., LTD.

Dressing accuracy of treatment
Dressing adhesive does not stick to hand
Dual removal technology of release paper
Secondary infection protective technology



OPILEX CO., LTD.